

# Laboratory organization system SmartRack®



e-Flyer



- 80 % Space Saving
- 10 % Increase in Efficiency
- Complies with workplace safety regulations in laboratories

Due to the current market situation with which we are all familiar, characterised by an exceptionally unstable pricing situation, the prices listed in the print medium may increase during the promotional period in individual cases. You can find the relevant applicable prices at [www.carlroth.com](http://www.carlroth.com)

**-15%**



A121°C

Made in Germany



Please note: Different conditions

## Laboratory organization system SmartRack®

Better Basics Laborbedarf. Material: aluminium powder-coated, aluminium anodized slotted rails, plastic (polyamide 12).

Modular organization system for any laboratory accessories and laboratory glassware. The wide array of different modules of the organization system can be set up in **any arrangement, combined** and used flexibly. Can adapt whenever and however the user requires. Plastic elements are temperature resistant, can be washed with mild cleaning agents or alcohol-based disinfectants. These are also resistant to chemicals and UV.

**Delivery incl.** frame, feet (1 pair), 3 slotted rails with side handles, collecting tray.

### SmartRack®

Width (mm)	Height (mm)	Depth (mm)	Art. No.	Pack Qty.	€	€
635	615	265	1PHE.1	1 unit	751,45	638,70



### Advantages:

- The modular design means users can adapt the system to suit their needs
- All modules can be combined in any quantity and any position
- 10% more efficiency thanks to one-handed operation and allows for optimised work processes
- Increases safety when working thanks to tilt protection
- 80% space savings thanks to improved use of workspace at height
- Creates a neater working environment with more free workspace

### SmartRack® Mini

Width (mm)	Height (mm)	Depth (mm)	Art. No.	Pack Qty.	€	€
420	380	265	1PHH.1	1 unit	536,45	455,95



**15%**

off the entire  
SmartRack® range

[www.carlroth.com](http://www.carlroth.com)





A 121 °C

Made in Germany

Please note: Different conditions

### Module set basic equipment for SmartRack®

Better Basics Laborbedarf. Material: aluminium powder-coated, plastic (polyamide 12).

Compatible with the SmartRack® laboratory organization system.



The basic module set is a collection of basic equipment for organising essential accessories for day-to-day laboratory use, and comprises 6 modules.

**Delivery incl.** 6 modules.

**Delivery without** SmartRack®, please order separately and without small accessories, waste disposal bags, tissue box, wash bottles.

- 1 x storage box (1PX6.1)
- 1 x storage box insert for small accessories (1PXK.1)
- 1 x holder for disposal bags (1PXL.1)
- 1 x holder for gloves/tissue boxes (1PXC.1)
- 1 x holder for wash bottles (1PXX.1)
- 1 x tripod adapter (1PX3.1)

Art. No.	Pack Qty.	€	€
1PHL.1	1 set	428,95	364,55



A 121 °C

Please note: Different conditions

### Holder for chromatography columns for SmartRack®

Better Basics Laborbedarf. Material: plastic (polyamide 12).

The module makes it easy to attach chromatography columns.

The holder adjusts to the column circumference and the snap-in mechanism provides a secure fit.

**Delivery without** chromatography columns.

Suitable for	Width (mm)	Art. No.	Pack Qty.	€	€
Chromatography column Ø 25–35 mm	75	1PX7.1	1 unit	75,25	63,90
Chromatography column Ø 36–46 mm	75	1PX8.1	1 unit	75,25	63,90



A 121 °C

Please note: Different conditions

### Holder for separating funnel for SmartRack®

Better Basics Laborbedarf. Material: aluminium powder-coated, plastic (polyamide 12).

The module makes it easy to attach separating/shaking funnels.

The Ø of the separating funnel is measured at its widest point.

**Delivery without** separating funnels.

Suitable for	Width (mm)	Art. No.	Pack Qty.	€	€
Ø 45–85 mm	75	1PTY.1	1 unit	64,50	54,75
Ø 85–105 mm	113	1PX0.1	1 unit	64,50	54,75
Ø ≥105 mm	133	1PX1.1	1 unit	64,50	54,75

**15 % off the entire SmartRack® range**

Carl Roth GmbH + Co. KG · Schoemperlenstr. 3-5 · 76185 Karlsruhe

Phone: +49 721 / 56 06 510 · Fax: +49 721 / 56 06 111 · E-Mail: info@carlroth.com · Internet: www.carlroth.com

For ordering and information about prices and delivery in your country please contact us or your local distributor.

Contact: www.carlroth.com/com/en/Auslandsservice · www.carlroth.com/com/en/Auslandsvertretungen

All supplies and deliveries are subject to the terms and conditions of sale and delivery of Carl Roth GmbH + Co. KG, Karlsruhe.