



Well advised with Roth.

Statement Amino Acid composition - typical data

1T92, Peptone ex soya **CELLPURE®** animal-free/GMO-free, for microbiology

The data below for the amino acid composition are only typical data and not part of the specification.

Amino Acids		Total (g/100g)
Alanine	Ala	3.18
Arginine	Arg	3.05
Aspartic acid	Asp	9.00
Cystine	Cys	0.79
Glutamic acid	Glu	11.83
Glycine	Gly	2.47
Histidine	His	1.23
Isoleucine	Ile	2.55
Leucine	Leu	3.90
Lysine	Lys	3.85
Metionine	Met	0.86
Phenylalanine	Phe	2.40
Proline	Pro	2.54
Serine	Ser	2.69
Treonine	Thr	2.61
Tryptophan	Trp	0.68
Tyrosine	Tyr	1.38
Valine	Val	2.94

Date of issue: 29th August 2022, Rev. 0

i.A. Jessica Wuske
Quality Management

