



Specification

Item number:	AE63		
Papanicolaou polychrome solution 3a (EA31)			
for microscopy			
CAS No.:		print date:	20.04.2017
formula:			
density:			
molcular weight:			

Type analysis

Suitability for microscopy complies

Appearance: muddy green-brown solution

Smell: methanole

Use test:

Staining with Papanicolaou-Polychromlösung EA 31 in dyeing step of Papanicolaou-staining.

Staining time:

- nuclear staining :	5 minutes
- Papanicolaou Orange II:	4 minutes
- Papanicolaou EA 31:	5 minutes

Staining result:

- Basophiles Zytoplasma	blue-green
- Eosinophiles Zytoplasma	rose
- Cell nuclei	blue

With Papanicolaou polychrome solution EA 31 gives a balanced cytoplasmic staining.

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

