

Specification

Item number: X906

Hematoxylin solution A acc. to Weigert for microscopy

CAS No.:

print date:

01.01.2024

formula:

density: 0,9

molecular weight:

Type analysis

Suitability for microscopy

complies

Appearance: red violet(Hämatoxylin)

Use test at histological sections:

For testing, prepare an additional solution Hematoxylin Weigert Solution B (ferric chloride solution).

Equally Solution A and B (according to the required amount) mix immediately before use.

Stain sections of distilled water 1 to 2 minutes .

May differentiate into hydrochloric acid-alcohol, rinse well in tap water.

Result:

Staining time: 1 - 2 minutes

Cell nuclei: blue-black

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

