

Spécification

Numéro d'article: X906

Hématoxyline, solution A selon Weigert

pour la microscopie

N° CAS:

Date d'impression: 01.01.2024

Formule:

Densité: 0,9

Masse moléculaire:

Type analysis

Suitability for microscopy complies

Appearance: red violet(Hämatoxylin)

Use test at histological sections:

For testing, prepare an additional solution Hematoxylin Weigert Solution B (ferric chloride solution).

Equally Solution A and B (according to the required amount) mix immediately before use.

Stain sections of distilled water 1 to 2 minutes .

May differentiate into hydrochloric acid-alcohol, rinse well in tap water.

Result:

Staining time: 1 - 2 minutes

Cell nuclei: blue-black

Nos produits sont contrôlés pour une utilisation en laboratoire.

Les indications se rapportent à l'état actuel de la qualité du produit.

Nous nous réservons le droit d'apporter des modifications si nécessaire.

Dr. R. Niemand
Head of Quality Assurance

S. Lorsee
Head of Quality Management

Ce document a été élaboré de manière informatique et ne comporte pas de signature.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

