





Valid for Item-No.: Tip 10 - 200µl = 739270, 739290, 739295
 Box with Tips = 739261, 739280




"All dimensions in mm"

Customer drawing subject to change without notice!

Uncontrolled Copy

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date	Date	Date	
-	01.03.2017	Mar. 2017	Mar. 2017	
Date	Name	Name	Name	
-	G. Kofler	Ch. Rathner	Ch. Rathner	

PDS No.739261/ 739280/ 739290	PRODUCT DATA SHEET			Page 1 of 1
Revision 03	Pipette Tips, 200 µl			
	Greiner Item-No. 739261 / 739280 / 739290			
Valid for Item-No.:	739261 (sterile)	739280	739290	

1.	Description / Specification	
1.1	Description	PP Pipette Tips, long tip collar 739261: sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	10 – 200 µl
1.4	Material / Resin	PP (Polypropylene), free of heavy metal Heavy metal-free colour
1.5	Colour	Yellow
1.6	Sterilization	739280, -290: no 739261: by irradiation
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Suited for Gilson® Pipetman, P-20, P-100, P-200

2.	Features	
2.1	Basic features	High breaking resistance (hard, stable form), autoclavable
2.2	Temperature range	- 190°C to approx. + 140°C
2.3	Autoclavability	121°C / 20 min.
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	739261: 3 years after month of production 739280, -290: n/a
2.7	Other Information	The testing procedures of the pipette tips follows the ISO 8655

3.	Packaging	739261, -280	739290
3.1	Pieces / Packaging unit	96 / yellow ST-rack (autoclavable polypropylene racks)	500 / Bulk pack (polyethylene bags)
3.2	Pieces / Box	4800	15000
3.3	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 22 April 2015	Date 23 April 2015	Date 23 April 2015	
Date 21.10.2013	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.