






Valid for Item-No.: 606190



"All dimensions in mm"  
Customer drawing subject to change without notice!

Uncontrolled Copy

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 1	Date 2011/05/11	Date 2016/05/18	Date 2016/05/18	
Date 2011/05/11	Name forb	Name forb	Name sonntag	

PDS No. 606190	<b>PRODUCT DATA SHEET</b>	Page 1 of 1
Revision 06	Pipette, CELLSTAR <sup>®</sup> , 5 ml, "shorty"	 <b>greiner bio-one</b>
	Greiner Item-No. 606190	

1.	Description / Specification	
1.1	Description	PS Pipette "shorty", graduation (scale 1/10 ml), blue nozzle, synthetic plug, 5 ml, sterile
1.2	Dimensions	See customer drawing
1.3	Volume	5 ml ; volumetric accuracy $\pm 2\%$
1.4	Material / Resin	Pipette and nozzle: PS (Polystyrene), free of heavy metal Plug: Bonded Polyester, free of heavy metal
1.5	Colour	Pipette: clear Nozzle: blue Plug: white
1.6	Sterilization	SAL $10^{-3}$
1.7	Quality Control	- Raw material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / bag	1 / "peel off"-packaging (plastic/plastic)
3.2	Pieces / case	200
3.3	Lot-No.	YYMMDDXX (year, month, day, index)
3.4	Other Information	Certificate of Quality

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 22 October 2014	Date 22 October 2014	Date 22 October 2014	
Date 15.07.2013	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.