



## TLC Plates

with Silica Gel 60 **-20%**



### Polyester foils POLYGRAM®

TLC-ready-to-use foil	Design	Plate format	Layer thickness	Art. No.	Pack Qty.	€	€
POLYGRAM® SIL G	-	4 x 8 cm	0,2 mm	N733.1	50	<del>34,90</del>	27,90
POLYGRAM® SIL G / UV <sub>254</sub>	with fluorescence indicator	4 x 8 cm	0,2 mm	0816.1	50	<del>34,50</del>	27,60
POLYGRAM® SIL G / UV <sub>254</sub>	with fluorescence indicator	4 x 8 cm	0,2 mm	<b>BIG PACK</b> 0816.2	200 (4 x 50)	<del>125,50</del>	100,40
POLYGRAM® SIL G / UV <sub>254</sub>	with fluorescence indicator	20 x 20 cm	0,2 mm	N737.1	25	<del>105,50</del>	84,40

### Aluminium foils ALUGRAM®/ALUGRAM® Xtra

TLC-ready-to-use foil	Design	Plate format	Layer thickness	Art. No.	Pack Qty.	€	€
ALUGRAM® SIL G	-	20 x 20 cm	0,2 mm	N728.1	25	<del>117,50</del>	94,00
ALUGRAM® SIL G / UV <sub>254</sub>	with fluorescence indicator	4 x 8 cm	0,2 mm	N729.1	50	<del>34,90</del>	27,90
ALUGRAM® SIL G / UV <sub>254</sub>	with fluorescence indicator	4 x 8 cm	0,2 mm	<b>BIG PACK</b> N729.2	200 (4 x 50)	<del>129,00</del>	103,20
ALUGRAM® SIL G / UV <sub>254</sub>	with fluorescence indicator	5 x 10 cm	0,2 mm	N730.1	50	<del>52,90</del>	42,30
ALUGRAM® SIL G / UV <sub>254</sub>	with fluorescence indicator	20 x 20 cm	0,2 mm	N732.1	20	<del>105,50</del>	84,40
ALUGRAM® Xtra SIL G / UV <sub>254</sub>	with fluorescence indicator	4 x 8 cm	0,2 mm	3874.1	50	<del>34,50</del>	27,60
ALUGRAM® Xtra SIL G / UV <sub>254</sub>	with fluorescence indicator	5 x 10 cm	0,2 mm	3876.1	50	<del>53,90</del>	42,90
ALUGRAM® Xtra SIL G / UV <sub>254</sub>	with fluorescence indicator	20 x 20 cm	0,2 mm	3879.1	25	<del>105,50</del>	84,40

Further DC-Plates and accessories in our current brochure Chemicals

# ROTI®CAP – GC Capillary Columns

MADE IN  GERMANY

## ROTI®CAP-5 MS Ultra Silarylene Phase

**-25%**

- **Extremely low column bleed**
- **Application areas:** all-round phase for environmental analyses, trace analyses, EPA methods, pesticides, PCBs, food and drug analyses. Ideal for ion-trap and quadrupol MS detectors
- **Similar phases:** DB-5 MS, HP-5 MS, Ultra-2, Equity™-5, CP-SIL 8 CB low bleed/MS, Rtx®-5SIL-MS, Rtx®-5 MS, Rxi®-5 MS, 007-5 MS, BPX™5, MDN-5S, AT™-5 MS, VF-5 MS
- Equivalent to USP phases G27/G36

Length (m)	Inner Ø (mm)	Film thickness (µm)	Art. No.	Pack Qty.	€	€
30	0,25	0,25	2959.1	1	<del>589,00</del>	441,75



**-20%**

## ROTI®CAP-5 MS

5 % Diphenyl - 95 % Dimethylpolysiloxane

- **Similar phases:** DB-5 MS, HP-5 MS, Ultra-2, Equity™-5, CP-SIL 8 CB low bleed/MS, Rtx®-5SIL-MS, Rtx®-5 MS, Rxi®-5 MS, 007-5 MS, BPX™5, MDN-5S, AT™-5 MS, VF-5 MS
- Equivalent to USP phases G27/G36

Length (m)	Inner Ø (mm)	Film thickness (µm)	Art. No.	Pack Qty.	€	€
30	0,25	0,25	CN15.1	1	<del>585,00</del>	468,00
30	0,32	0,25	1478.1	1	<del>629,00</del>	502,90

## ROTI®CAP-5

5 % Phenyl - 95 % Methylpolysiloxane

- **Similar phases:** SE-54, SE-52, DB-5, HP-5, SPB™-5, CP-SIL 8, Rtx®-5, 007-5, BP5, MDN-5, AT™-5, ZB-5
- Equivalent to USP phases G27/G36

Length (m)	Inner Ø (mm)	Film thickness (µm)	Art. No.	Pack Qty.	€	€
30	0,25	0,25	CN11.1	1	<del>500,00</del>	406,90
30	0,25	1,00	1386.1	1	<del>500,00</del>	406,90
30	0,32	0,25	1227.1	1	<del>539,00</del>	430,90
30	0,32	0,50	1389.1	1	<del>539,00</del>	430,90

## ROTI®CAP-624

6 % Cyanopropyl-phenyl - 94 % Dimethylpolysiloxane

- **Similar phases:** HP-624, HP-VOC, DB-624, DB-VRX, SBP™-624, CP-624, Rtx®-624, Rtx®-Volatiles, 007-624, BP624, VOCOL®, ZB-624
- Equivalent to USP phase G43

Length (m)	Inner Ø (mm)	Film thickness (µm)	Art. No.	Pack Qty.	€	€
30	0,32	1,80	CN16.1	1	<del>540,00</del>	438,90
60	0,32	1,80	1876.1	1	<del>969,00</del>	775,20

## ROTI®CAP-WAX

Polyethylene glycol (PEG) 20000 dalton

- **Similar phases:** DB-Wax, Supelcowax®, HP-Wax, HP-INNOWAX, Rtx®-Wax, CP-Wax 52 CB, Stabilwax®, 007-CW, BP20, AT™-Wax, ZB-Wax
- Equivalent to USP phase G16

Length (m)	Inner Ø (mm)	Film thickness (µm)	Art. No.	Pack Qty.	€	€
25	0,53	1,00	6386.1	1	<del>629,00</del>	502,90
30	0,25	0,25	CN12.1	1	<del>525,00</del>	420,00
30	0,32	0,25	1618.1	1	<del>539,00</del>	430,90
30	0,32	0,50	1228.1	1	<del>539,00</del>	430,90

## ROTI®CAP-FFAP

Polyethylene glycol 2-nitroterephthalate

- **Similar phases:** DB-FFAP, HP-FFAP, CP-SIL 58 CB, 007-FFAP, CP-FFAP CB, Nukol™
- Equivalent to USP phase G35, close equivalent to G25

Length (m)	Inner Ø (mm)	Film thickness (µm)	Art. No.	Pack Qty.	€	€
30	0,25	0,25	CN13.1	1	<del>500,00</del>	406,90
30	0,32	0,25	1229.1	1	<del>550,00</del>	446,90
30	0,32	0,50	1854.1	1	<del>539,00</del>	430,90

The desired column is not present?

For our extended assortment of **capillary columns**, see chapter **Chromatography** in our current catalogue or at [www.carlroth.com](http://www.carlroth.com)

We will gladly submit a quote!

### Contact international

Carl Roth GmbH + Co. KG  
Tel: +49 (0)721 5606-0  
Fax: +49 (0)721 5606-149

Schoemperlenstr. 3-5  
D-76185 Karlsruhe  
[www.carlroth.com](http://www.carlroth.com)

### Contact Netherlands

Tel.: 0180/51 67 04  
Fax: +49 721/56 06-260  
[bestellingen@carlroth.nl](mailto:bestellingen@carlroth.nl)  
[www.carlroth.nl](http://www.carlroth.nl)

Carl Roth GmbH + Co. KG · Schoemperlenstr. 3-5 · D-76185 Karlsruhe  
De algemene verkoop- en leveringsvoorwaarden van Carl Roth GmbH + Co. KG, Karlsruhe zijn van toepassing.  
Les conditions générales de vente et de livraison de Carl Roth GmbH + Co. KG, Karlsruhe sont en vigueur.

### Contact Belgium

Tel.: 03 283 47 10 (Vlaanderen & Brussel) · Fax: +49 721/5606-260  
Tél.: 080 447 958 (Wallonie & Bruxelles) · Fax: +49 721/5606-259  
[bestellingen@carlroth.be](mailto:bestellingen@carlroth.be) · [commandes@carlroth.be](mailto:commandes@carlroth.be)  
[www.carlroth.be](http://www.carlroth.be)