

Phosphate buffered saline with or without Tween® for cell biology, molecular biology and biochemistry.

## ... as solution



| Product name        | Purity  | Endotoxin content | pH value  | Art. No. | Pack Qty. |
|---------------------|---|-------------------|-----------|----------|-----------|
| ROTI®Cell DPBS      | ready-to-use, sterile, w/o Ca/Mg                            | ≤1,0 EU/ml        | 7,3 ±0,3  | 9124.1   | 500 ml    |
|                     |   |                   |           | 9124.2   | 1 l       |
| ROTI®Cell 10x DPBS  | 10x conc., sterile, w/o Ca/Mg                               | ≤1,0 EU/ml        | 6,8 ±0,2  | 9130.1   | 500 ml    |
|                     |   |                   |           | 9130.2   | 1 l       |
| ROTI®Cell DPBS      | ready-to-use, sterile, with Ca/Mg                           | ≤1,0 EU/ml        | 7,0 ±0,2  | 9131.1   | 500 ml    |
| ROTI®Cell PBS       | ready-to-use, sterile, w/o Ca/Mg                            | ≤1,0 EU/ml        | 7,4 ±0,1  | 9143.1   | 500 ml    |
|                     |   |                   |           | 9143.2   | 1 l       |
| ROTI®Cell 10x PBS   | 10x conc., sterile, w/o Ca/Mg                               | ≤1,0 EU/ml        | 7,4 ±0,1  | 9150.1   | 1 l       |
| ROTI®Stock 10x PBS  | 10x conc., BioScience Grade, ready-to-use, sterile filtered |                   | 7,4 ±0,25 | 1058.1   | 1 l       |
| ROTI®Stock 10x PBST | 10x conc., BioScience Grade, ready-to-use, sterile filtered |                   | 7,4 ±0,25 | 1059.1   | 1 l       |

For safety information and additional data, see our current catalogue or at [www.carlroth.com](http://www.carlroth.com)

## ... as tablet



| Product name                     | Batch quantity | pH value  | Art. No.  | Pack Qty.   |             |
|----------------------------------|----------------|-----------|-----------|-------------|-------------|
| ROTI®Fair PBS 7.2                | 1000 ml/tablet | 7,2 ±0,05 | 1106.1    | 10 unit(s)  |             |
|                                  |                |           | 1106.2    | 100 unit(s) |             |
| ROTI®Fair PBS 7.4                | 100 ml/tablet  | 7,4 ±0,05 | 1107.1    | 100 unit(s) |             |
|                                  | 200 ml/tablet  | 7,4 ±0,05 | 1108.1    | 100 unit(s) |             |
|                                  | 500 ml/tablet  |           | 7,4 ±0,05 | 1111.1      | 12 unit(s)  |
|                                  |                |           |           | 1111.2      | 100 unit(s) |
|                                  | 1000 ml/tablet |           | 7,4 ±0,05 | 1112.1      | 10 unit(s)  |
|                                  |                |           |           | 1112.2      | 100 unit(s) |
| ROTI®Fair PBS potassium-free 7.4 | 1000 ml/tablet | 7,4 ±0,05 | 1113.1    | 10 unit(s)  |             |
|                                  | 1000 ml/tablet | 7,4 ±0,05 | 1114.1    | 100 unit(s) |             |
| ROTI®Fair PBST 7.4               | 500 ml/tablet  | 7,4 ±0,05 | 1115.1    | 12 unit(s)  |             |
|                                  |                |           | 1115.2    | 100 unit(s) |             |
|                                  | 1000 ml/tablet | 7,4 ±0,05 | 1116.1    | 10 unit(s)  |             |
|                                  |                |           | 1116.2    | 100 unit(s) |             |

For safety information and additional data, see our current catalogue or at [www.carlroth.com](http://www.carlroth.com)



## ... or as powder

| Product name          | Batch quantity | pH value  | Art. No. | Pack Qty. |
|-----------------------|----------------|-----------|----------|-----------|
| ROTI®Fair PBS 7.4     | 10 l/pouch     | 7,4 ±0,05 | 1101.1   | 1 unit(s) |
|                       | 50 l/pouch     | 7,4 ±0,05 | 1103.1   | 1 unit(s) |
|                       | 100 l/pouch    | 7,4 ±0,05 | 1104.1   | 1 unit(s) |
| ROTI®Fair 10x PBS 7.4 | 1000 ml/pouch  | 7,4 ±0,05 | 1105.1   | 5 unit(s) |
| ROTI®PreMix PBS       | 9.55 g/l       | 7,4 ±0,2  | 0890.1   | 250 g     |
|                       |                |           | 0890.2   | 1.25 kg   |
| ROTI®PreMix PBST      | 10.05 g/l      | 7,4 ±0,2  | 0987.1   | 250 g     |
|                       |                |           | 0987.2   | 1.25 kg   |

For safety information and additional data, see our current catalogue or at [www.carlroth.com](http://www.carlroth.com)

Current prices at [www.carlroth.com](http://www.carlroth.com)

Contact international:

Phone +49 721/5606 510 · Fax +49 721/5606 111 · [info@carlroth.com](mailto:info@carlroth.com) · [www.carlroth.com](http://www.carlroth.com)

Carl Roth GmbH + Co. KG · Schoemperlenstr. 3-5 · D-76185 Karlsruhe

Export prices might be higher.

All supplies and deliveries are subject to the terms and conditions of sale and delivery of Carl Roth GmbH + Co. KG, Karlsruhe