



ROTH-Media

acc. to European Pharmacopoeia

Issue 6.0 (2008), 6.3 (2009), 7.0 (2011) and 8.0 (2014)

Medium	Art. No.	acc. to Pharm. Eur. 6.0 (2008)	acc. to Ph. Eur. 6.0 (2008, Harmonised Method), 6.3 (2009), 7.0 (2011) or 8.0 (2014)
Antibiotic Medium No. 1	CL53	<i>not listed in this norm</i>	Antibiotic Agar No. 1 (Seed Agar) <i>(not listed in 8.0)</i>
Baird Parker Agar (Base)	X913	Agar Medium O	<i>not listed in this norm</i>
Brilliant Green Tetrathionate Bile Broth	AE42	Broth Medium I	<i>not listed in this norm</i>
Brilliant Green Lactose Sucrose Agar	AE43	Agar Medium L	<i>not listed in this norm</i>
CASO Agar	X937	Agar Medium B	Casein Soya Bean Digest Agar
CASO Broth	X938	Broth Medium A	Casein Soya Bean Digest Broth
CASO Broth, granulated	6678	Broth Medium A	Casein Soya Bean Digest Broth
Cetrimide Agar (Base)	X918	Agar Medium N	Cetrimide Agar
Reinforced Clostridial Medium	AE66	Medium P	Reinforced Medium for <i>Clostridia</i>
Columbia Agar (Base)	X919	Medium Q	Columbia Agar
Desoxycholate Citrate Agar	AE44	Agar Medium J	<i>not listed in this norm</i>
Potato Extract Glucose Agar	X931	<i>not listed in this norm</i>	Potato Dextrose Agar
Violet Red Bile Agar with Glucose and Lactose	AE45	Agar Medium F	<i>not listed in this norm</i>
Violet Red Bile Agar with Glucose	X939	<i>not listed in this norm</i>	Violet Red Bile Glucose Agar
Lactose-Broth	X921	Broth Medium D	<i>not listed in this norm</i>
Lactose Sulfite Broth (Base)	AE46	Medium R	<i>not listed in this norm</i>
MacConkey Agar	X922	Agar Medium H	MacConkey Agar
MacConkey Broth	AE47	Broth Medium G	MacConkey Broth
Mannitol Salt Agar	CL81	<i>not listed in this norm</i>	Mannitol Salt Agar
Mossel EE Broth	AE65	Enriched Broth Medium E	Enterobacteria enrichment broth-Mossel
Peptone Water, buffered	AE67	Buffered Sodium Chloride-Peptone Solution pH 7.0	Buffered Sodium Chloride-Peptone Solution pH 7.0
R2A Agar	CL01	Agar Medium S	<i>not listed in this norm</i>
Rappaport Vassiliadis Broth	CL03	<i>not listed in this norm</i>	Rappaport Vassiliadis <i>Salmonella</i> Enrichment Broth
Sabouraud 2 % Glucose Bouillon	AE23	<i>not listed in this norm</i>	Sabouraud Dextrose Broth
Sabouraud 4% Glucose Agar	X932	Agar Medium C	Sabouraud Dextrose Agar
Sabouraud 4% Glucose Agar (with Chloramphenicol)	AE48	Agar Medium C	<i>not listed in this norm</i>
Triple Sugar Iron Agar	AE49	Agar Medium M	<i>not listed in this norm</i>
XLD Agar	X941	Agar Medium K	Xylose, Lysine, Desoxycholate Agar

s.s. 12/2015