



Specification

Item number:	5783		
Roti®-Load DNASTain 1			
6x conc. gel loading buffer, ready-to-use, DNase-free			
CAS No.:		print date:	22.04.2017
formula:			
density:		storage temperature:	-20 °C
molcular weight:			

Appearance: dark blue, clear solution
Suitable for DNA fragments: > 500 bp

Typical formulation:

Tris-HCl (pH 7,5): 60 mM
EDTA: 60 mM
Glycerol: 60% (w/v)
EvaGreen®DNA stain
Bromphenol blue
Xylene cyanole

All components of the buffer are proved to be DNase-free.
The production is carried out under sterile conditions.

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

