



# Specification

---

Item number: 0124

## PTFE paste

CAS No.:

print date:

01.01.2024

formula:

density:

molecular weight:

---

Suitable for sealing fine threads and cuts. High viscosity, good lubrication properties, good electrical resistance properties. Sealed thread or quick fit lid can be opened easily. Good resistance against acids, bases and standard solvents - except fluor components. Temperature stable from -30 to +290 °C.

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

---

**Dr. R. Niemand**

*Head of Quality Assurance*

**S. Lorsee**

*Head of Quality Management*

*This document is computer printout and has therefore not been signed by hand.*

---

**Carl Roth GmbH + Co. KG**

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: [info@carlroth.de](mailto:info@carlroth.de)

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

