



# Specification

---

Item number: 0984

## ROTI®GelStain Red

**20 000x conc., ready-to-use**

CAS No.:

print date:

01.01.2026

formula:

density:

storage temperature: +4 °C

molecular weight:

---

### *Type analysis*

#### **Appearance**

red, transparent liquid, no precipitation

#### **Staining of DNA in agarose gels**

complies

#### **Staining of DNA in PA gels**

complies

Excitation maximum (bound to DNA)

approx. 300 nm and 530-550 nm

Emission maximum (bound to DNA)

630 nm

Each batch is tested for display in the electrophoresis.

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

---

**Dr. R. Niemand**

*Head of Quality Assurance*

**S. Lorsee**

*Head of Quality Management*

*This document is computer printout and has therefore not been signed by hand.*

---

**Carl Roth GmbH + Co. KG**

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: [info@carlroth.de](mailto:info@carlroth.de)

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet