



Specification

Item number: 1057

ROTI®Stock 20 % SDS

20 % SDS, BioScience Grade, ready-to-use, sterile filtered

CAS No.:

print date:

01.01.2024

formula:

density: ~1

storage temperature: +15 bis +25 °C

molecular weight:

Type analysis

Appearance of solution

colorless liquid

Colour

water-clear

20 % SDS in deionised water, 0.2 µm filtered

Store the solution away from light in a dry place. SDS stock solution and other solutions containing SDS can precipitate easily at room temperature. The crystals can be easily resolubilised by brief heating them between 30 to 50 °C. Prolonged heating above 30 °C should be avoided, otherwise hydrolysis of the SDS and formation of fatty alcohols may occur.

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

