

Specification

Item number: 1L8H

Peptone Water, buffered

ISO 6579, ISO 6887, ISO 22964, for microbiology

CAS No.:

print date:

01.01.2022

formula:

density: 1,03

storage temperature: +15 bis +25 °C

molecular weight:

Casein peptone	10 g/l
Sodium chloride (NaCl)	5 g/l
Disodium hydrogen phosphate	3.5 g/l
Dipotassium hydrogen	1.5 g/l
pH-value (25 °C)	6.8-7.2

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

