

Specification

Item number: 1NN1

anti-alpha-Smooth Muscle Actin mouse monoclonal, 1A4/ASM-1

100 µg/ml, monoclonal, affinity purified

CAS No.:

print date:

01.01.2025

formula:

density: ~1

storage temperature: +4 °C

molecular weight:

Type analysis

Clonality	Monoclonal
Concentration	100 µg/ml
Raised in	Hybridoma cell culture supernatant
Immunglobulin (Ig) class	IgG2a
Clone	1A4/ASM-1
Conjugate	Unconjugated
Purification	Affinity chromatography
Purity (SDS-PAGE)	≥95 %

Immunogen: Synthetic N-terminus decapeptide of alpha-smooth-muscle isoform of actin

PBS buffer, pH 7.4 with 0.09 % sodium azide and 0.5 % BSA.

Short term at 2-8 °C (+4 °C);

long term storage in aliquots at -20 °C; avoid freeze/ thaw cycles.

Manufacturer: PROGEN

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

