

Specification

Item number: 1NN6

anti-Desmocollin 1 mouse monoclonal, Dsc1-U100

50 µg/ml, monoclonal, affinity purified

CAS No.:

print date:

01.01.2024

formula:

density: ~1

storage temperature: +4 °C

molecular weight:

Purification	Affinity chromatography
Clonality	Monoclonal
Concentration	50 µg/ml
Raised in	Hybridoma cell culture supernatant
Immunglobulin (Ig) class	IgG1
Clone	Dsc1-U100
Conjugate	Unconjugated
Purity (SDS-PAGE)	≥95 %

Immunogen: Synthetic peptide corresponding to sequence present within the intracellular part of human desmocollin 1 (HQTLESVKGVGGDTGRYATDWQC)

PBS buffer, pH 7.4 with 0.09 % sodium azide and 0.5 % BSA.

Short term at 2-8 °C (+4 °C);

long term storage in aliquots at -20 °C; avoid freeze/ thaw cycles.

Manufacturer: PROGEN

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

