

Specification

Item number: 1NN7

anti-HA-tag mouse monoclonal, 12CA5

lyophilized, monoclonal, affinity purified

CAS No.:

print date:

01.01.2024

formula:

density:

storage temperature: +4 °C

molecular weight:

Purification	Affinity chromatography
Clonality	Monoclonal
Raised in	Hybridoma cell culture supernatant
Immunglobulin (Ig) class	IgG2b
Clone	12CA5
Conjugate	Unconjugated
Purity (SDS-PAGE)	≥95 %

Immunogen: 9 amino acids from the human influenza virus hemagglutinin protein

Lyophilized; reconstitute in sterile PBS, pH 7.4.

Storage before reconstitution at 2-8 °C until indicated expiry date.

Storage after reconstitution at -20 °C (avoid freeze/thaw cycles).

Manufacturer: PROGEN

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

