

Specification

Item number: 2493

ROTI®DipSlide PCA/Legionella

ISO 11133, ready-to-use, sterile, for microbiology

CAS No.:

print date:

01.01.2024

formula:

density:

storage temperature: +4 bis +15 °C

molecular weight:

Composition

Agar side 1:

Casein peptone	5 g/l
Yeast extract	2,5 g/l
Glucose	1 g/l
Agar	15 g/l
pH-value	7,0±0,2
Colour	amber, slightly opalescent
Growth of test strains	complies

Agar side 2:

Yeast extract	10 g/l
Charcoal activated	2 g/l
GVPC selective additive	complies
Growth additive	complies
Agar	13 g/l
pH-value	6,9±0,2
Colour	black
Growth of test strains	complies

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

