

Specification

Item number: 3129

ROTI®Fect RNAi Kit

CELLPURE® ready-to-use, for transfection

CAS No.:

print date:

01.01.2024

formula:

density:

storage temperature: -20 °C

molecular weight:

Roti®-Fect RNAi Lipo and Roti®-Fect RNAi Puffer (10x conc.)

Appearance Roti®-Fect RNAi Lipo clear solution

Appearance Roti®-Fect RNAi Puffer clear solution

Ready-to-use, aqueous liposoma formulation and 10x Puffer, special developed for transfection of mammalian cells with siRNA or miRNA.

Successfully tested for following cell types:

Primary cells: MCS, M1

Cell lines: A549, 3T3-L1, B16F10, CHO-K1, H441, HT1080, HeLa, HeLa-KB, HepG2

For approx.250 transfections per ml im 24 wells.

Sterility test for Thioglycolate medium (7d / 37 °C)

Each batch test of transfection efficiency.

Not suitable for diagnostic, therapeutic or other clinical use on humans or animals.

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

