



Specification

Item number: 3735

Anti His Serum

polyclonal

CAS No.:

print date: 21.04.2017

formula:

density:

storage temperature: -20 °C

molcular weight:

Type analysis

Appearance	clear, ochre-rose coloured liquid
Raised in	rabbit
Clonality	polyclonal
Target	His-tag
Western Blot	complies
ELISA	complies

Recommended Application: Immunoblotting (Western-, Dotblot), ELISA, Proteinisolation

Purification: nothing

Buffer Solution: nothing

Stabilisers: nothing

Each batch is tested for Western-Blot and ELISA.

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

