



Spécification

Numéro d'article: 3735

Anti-His sérum

polyclonal

N°CAS

Date d'impression: 21.04.2017

Formule:

Densité:

Temp. de stockage: -20 °C

Masse moléculaire:

Type analysis

Appearance	clear, ochre-rose coloured liquid
Raised in	rabbit
Clonality	polyclonal
Target	His-tag
Western Blot	complies
ELISA	complies

Recommended Application: Immunoblotting (Western-, Dotblot), ELISA, Proteinisolation

Purification: nothing

Buffer Solution: nothing

Stabilisers: nothing

Each batch is tested for Western-Blot and ELISA.

Nos produits sont contrôlés pour une utilisation en laboratoire.
Les indications se rapportent à l'état actuel de la qualité du produit.
Nous nous réservons le droit d'apporter des modifications si nécessaire.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

Ce document a été élaboré de manière informatique et ne comporte pas de signature.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carloth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

