

Specification

Item number:	6933		
Boc-L-Serine			
≥99 %, for biochemistry			
CAS No.:	3262-72-4	print date:	10.11.2020
formula:	C ₈ H ₁₅ NO ₅		
density:		storage temperature:	+4 °C
molecular weight:	205,21 g/mol		

Type analysis

Appearance	white to off-white powder
Purity (HPLC)	≥99,0 %
Assay (NaOH)	98,0-102,0 %
Specific rotation $[\alpha]_D^{20}$ (1 % in H₂O)	-9,5° bis -7,5°
Water (KF)	≤0,5 %
Enantiomer purity	≤0,5 % D-Enantiomer

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

Dr. R. Niemand
Head of Quality Assurance

S. Lorsee
Head of Quality Management

This document is computer printout and has therefore not been signed by hand.