

BCMI-004 | CuSO4 (baseclick grade)

For research only!

Baseclick No: BCMI-004

Product name: CuSO₄ (baseclick grade)

CAS No: 7758-98-7

Molecular formula: CuSO₄

Molecular weight: 159.6 g/mol

Purity: ≥ 98 % (iodometric)

Appearance: light green to grey fine powder

Solubility: aqueous media

Storage Conditions: RT, dry, inert gas

Shelf life: 12 months unopened after receipt

Application

Water-soluble copper(I) source for click reactions after reduction (e.g. by using BCMI-005)

For References and FAQs see Online