

BCAA-001 | 4-Azido-L-phenylalanine

For research only!

Baseclick No: BCAA-001

Product name: 4-Azido-L-phenylalanine

CAS No: 33173-53-4

Molecular formula: $C_9H_{10}N_4O_2$

Molecular weight: 206.2 g/mol

Purity: ≥ 98 % (HPLC)

Appearance: off-white to brownish powder

Solubility: c = 5% in AcOH 80% (clear slightly colored solution)

Storage Conditions: -20 °C, dry

Shelf life: 12 months unopened after receipt

 $H_2N_{1,1}$ OH

For References and FAQs see Online