

Specification

Item number: 8272

200 bp-DNA-Ladder

ready-to-use, not prestained

CAS No.:

print date:

01.01.2024

formula:

density: 1

storage temperature: -20 °C

molecular weight:

Appearance blue, clear solution

Typical Composition:

TRIS/HCl (pH 7.5) 10 mM

EDTA 10 mM

Glycerol 10 % (w/v)

Bromphenol blue

Xylene cyanole

DNA Concentration 105 ng/μl

Size Range 200 to 2 kbp

Numbers of DNA Fragments 10

Ready-made solution

Generated by PCR amplification based on lambda phage DNA.

Fragment size at 200, 400, 600, 800, 1000, 1200, 1400, 1600, 1800 and 2000 bp.

Reinforced band at 1000 bp.

Each batch is tested for functionality in the electrophoresis.

Instructions for use are attached.

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

