

# Specification

Item number: 9248

## ROTI®Pol Hot-TaqS Mix (2x)

**2x conc., ready-to-use**

CAS No.:

print date:

01.01.2022

formula:

density: 1

storage temperature: -20 °C

molecular weight:

### Type analysis

<b>Appearance</b>	colourless
<b>Endonuclease activity</b>	none detected
<b>Exonucleases activity</b>	none detected
<b>Suitability for PCR (400 bp)</b>	complies
<b>Suitability for PCR (3 kb)</b>	complies
<b>Suitability for PCR (5 kb)</b>	complies
<b>Selective specificity in PCR</b>	complies

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

**Dr. R. Niemand**

*Head of Quality Assurance*

**S. Lorse**

*Head of Quality Management*

*This document is computer printout and has therefore not been signed by hand.*

**Carl Roth GmbH + Co. KG**  
Schoemperlenstraße 3-5  
76185 Karlsruhe

Telefon 0721/5606-0  
Telefax 0721/5606-149  
E-Mail: [info@carlroth.de](mailto:info@carlroth.de)

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

