

Specification

Item number: 9650

ROTI®Aquatest Plate Staph

Ph. Eur., ready-to-use, sterile, for microbiology

CAS No.:

print date: 01.01.2024

formula:

density:

storage temperature: +4 bis +15 °C

molecular weight:

Composition

Pancreatic digest of Casein	5 g/l
Peptone (peptic digest of animal tissue)	5 g/l
Beef extract	1 g/l
Sodium chloride (NaCl)	75 g/l
D-Mannitol	10 g/l
Phenol red	0.025 g/l
Agar	15 g/l
pH-value	7.4 ±0.2
Colour	red
Growth of test strains	complies
Diameter	55 mm

Enumeration of staphylococci.

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

