

# Specification

---

Item number: P078

**ROTI®Lumin**

**10x conc., fluorescence-free, for protein and cytoimmunochemistry**

CAS No.:

print date:

01.01.2024

formula:

density:

storage temperature: +4 °C

molecular weight:

---

A kit consists of:

120 ml (25 ml) Roti®-Lumin 1 and 120 ml (25 ml) Roti®-Lumin 2.

Sufficient for approx. 2400 cm<sup>2</sup> (500 cm<sup>2</sup>) membrane.

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

---

**Dr. R. Niemand**

*Head of Quality Assurance*

**S. Lorsee**

*Head of Quality Management*

*This document is computer printout and has therefore not been signed by hand.*

---

**Carl Roth GmbH + Co. KG**

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: [info@carlroth.de](mailto:info@carlroth.de)

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

