



Specification

Item number:	P730		
Roti®-Store cryo vials			
CAS No.:		print date:	02.07.2020
formula:			
density:		storage temperature:	+4 °C
molcular weight:			

Appereance: clear, yellow solution with 25 blue glass beads
Size: Length 47 mm, Diameter 11 mm
Sterilization at Ph.Eur. at least 15 min. at 121°C
Sterilization test no growth at 30°C after 3 days/n o growth at 20°C after 2 days
pH-value 7,15 - 7,65

Ready-to-use

Medium for the cryopreservation of microorganisms in general

Cryomedium composition: Rich growth medium for microorganisms and anti-freeze additive

Cryo long-term storage: -80 °C or liquid N2

Cryo short term storage: -20 °C

Successfully tested on clostridia, corynebacteria, lactobacilli, Staphylococci, Enterobacteria

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

